

Objective: Study the materials required to stop radiation from being transmitted.

State Standard: SCI.CP.7.7 2010

**Setup**

- 1) Turn the ST-160 nuclear lab station on, the button is on the back of the unit.
- 2) Set the voltage on the Geiger tube by pressing the H.V. button and push the up button until the display reads 400.
- 3) Push the H.V. button again; the red LED should go out.

**Measure the background radiation level**

- 1) Obtain a timer from your teacher.
- 2) Push the count button and start the timer, let the counter run for 30 seconds.
- 3) Push stop after 30 seconds.
- 4) Record the counts and repeat 3 times and take the average.
  - a. Run 1 \_\_\_\_\_
  - b. Run 2 \_\_\_\_\_
  - c. Run 3 \_\_\_\_\_
  - d. Average \_\_\_\_\_

**Testing Alpha particles**

- 1) Get a set of sources from your teacher. These samples are radioactive but at a very low level and are safe to handle.
- 2) Keep the sources at least 2 feet from the ST-160.
- 3) Insert the alpha source (Po-210), label down, into the plastic sample tray and place into the 3<sup>rd</sup> slot from the top inside the test opening.
- 4) Conduct three sets of counts of 30 seconds each and average the results.
  - a. Run 1 \_\_\_\_\_
  - b. Run 2 \_\_\_\_\_
  - c. Run 3 \_\_\_\_\_
  - d. Average \_\_\_\_\_
- 5) Insert the 4 mil poly shield (#1) into the top slot and conduct 3 counts and average as before.
  - a. Run 1 \_\_\_\_\_
  - b. Run 2 \_\_\_\_\_
  - c. Run 3 \_\_\_\_\_
  - d. Average \_\_\_\_\_
- 6) How much radiation is reaching the detector through the shield?
  - a.  $\frac{\text{Average from step 5}}{\text{Average from step 4}} \times 100 = \text{_____} \%$

- 7) Remove the Poly shield (#1) and replace it with the .04" plastic shield (#4). Conduct another 3 counts and average.
  - a. Run 1 \_\_\_\_\_
  - b. Run 2 \_\_\_\_\_
  - c. Run 3 \_\_\_\_\_
  - d. Average \_\_\_\_\_
- 8) How much radiation is reaching the detector through the shield?
  - a.  $\frac{\text{Average from step 7}}{\text{Average from step 4}} \times 100 = \text{_____}\%$
- 9) Remove the plastic shield and source and return them to their storage locations.

**Testing Beta Particles**

1. Insert the beta source (Sr-90), label down, into the plastic sample tray and place into the 3<sup>rd</sup> slot from the top inside the test opening.
2. Conduct three sets of counts of 30 seconds each and average the results.
  - a. Run 1 \_\_\_\_\_
  - b. Run 2 \_\_\_\_\_
  - c. Run 3 \_\_\_\_\_
  - d. Average \_\_\_\_\_
3. Insert the .04 plastic shield (#4) into the top slot and conduct 3 counts and average as before.
  - a. Run 1 \_\_\_\_\_
  - b. Run 2 \_\_\_\_\_
  - c. Run 3 \_\_\_\_\_
  - d. Average \_\_\_\_\_
4. How much radiation is reaching the detector through the shield?
  - a.  $\frac{\text{Average from step 3}}{\text{Average from step 2}} \times 100 = \text{_____}\%$
5. Remove the .04 plastic shield (#4) and replace it with the .09" aluminum shield (#7). Conduct another 3 counts and average.
  - a. Run 1 \_\_\_\_\_
  - b. Run 2 \_\_\_\_\_
  - c. Run 3 \_\_\_\_\_
  - d. Average \_\_\_\_\_
6. How much radiation is reaching the detector through the shield?
  - a.  $\frac{\text{Average from step 5}}{\text{Average from step 2}} \times 100 = \text{_____}\%$
7. Remove the .09" aluminum shield (#7) and replace it with the .032" lead shield (#8). Do three counts and average as before.
  - e. Run 1 \_\_\_\_\_
  - f. Run 2 \_\_\_\_\_
  - g. Run 3 \_\_\_\_\_
  - h. Average \_\_\_\_\_
8. How much radiation is reaching the detector through the shield?



**Analysis**

1) What particle was the easiest to stop?

---

---

2) Which source was the most active? How can you tell?

---

---

---

---

3) Which source was the hardest to stop?

---

---

---

---

4) How would you protect yourself if you had to go into an environment that was contaminated with alpha emitting radiation?

---

---

---

---

---